



# Patient Positioning in the Intensive Care Unit

*HoverMatt • HoverSling • Q2Roller • Positioning Wedges • HoverCover • HoverJack*

# ELEVATE YOUR STANDARD OF CARE

Elevating staff safety, patient care, and value



## ELEVATING STAFF SAFETY

**PROBLEM:** On average, nurses move the weight of a rhinoceros per 8-hour shift.<sup>1</sup> Hospital staff are injured at 2x the rate of private industry, with 5.7 injuries per 100 FTE's.<sup>2</sup>

**SOLUTION:** HoverMatt reduces the force needed to move a patient by 80–90%.<sup>3,4</sup> Evidence shows hospitals **reduced patient-handling-related staff injuries by 73%, and lost and restricted workdays by 94%**<sup>5</sup> by using the HoverMatt in conjunction with an effective Safe Patient Handling and Mobility Program.



## ELEVATING PATIENT CARE

**PROBLEM:** 60,000 patients die due to pressure injuries each year.<sup>6</sup> About 1 in 10 patients incur pressure injuries, resulting in more care and longer inpatient stays.<sup>7</sup>

**SOLUTION:** Supporting your efforts to decrease pressure injury rates, the HoverMatt **does not impair the microclimate or effect localized pressure, nor increase the risk to skin and soft tissue due to friction and shear.**<sup>8</sup>



## ELEVATING VALUE

**PROBLEM:** The average cost to replace a nurse is \$27K–\$103K.<sup>9</sup> A workers compensation claim from patient-handling-related injuries averages \$56,400 in direct and indirect costs.<sup>10,11</sup> With Nursing shortages expected to intensify and insufficient staffing, nurse stress levels will likely increase as job satisfaction and patient care decrease.<sup>12</sup>


**SOLUTION:** Evidence shows hospitals **reduced direct medical and indirect costs by 98%** by using the HoverMatt in conjunction with an effective Safe Patient Handling and Mobility Program.<sup>5</sup> **50% or less staff are needed to safely move a patient with HoverMatt products.**<sup>13</sup>



# BREATHABILITY



**Breathability:** The ability of a fabric to allow moisture vapor to be transmitted through the material, AKA: the Moisture Vapor Transmission Rate (MVTR). The higher the MVTR, the more breathable the fabric.



The **HoverMatt® Single-Patient Use Air Transfer Mattress** exceeds the breathability threshold by 23%.

- *HoverMatt's construction allows it to safely remain under the patient for length of stay.*



The **HoverSling® Repositioning Sheet** fabric pulls heat and moisture away from the patient, creating the ideal microclimate.

- *Third party testing proves this device does not affect the efficacy of a low air loss mattress, giving you the confidence to leave it under your patients for all-day care.*

**HoverCover™  
Breathable Absorbent  
Pad**

**Pro Tip:** Pair with any HoverMatt or HoverSling product for exceptional fluid protection and easy cleanup.



# HOVERMATT® SPU

Transfer and Position with Safety and Ease

Designed to be multifunctional, the breathable HoverMatt® Single-Patient Use (SPU) Air Transfer System makes patient transfers, turning, boosting, positioning and proning easier, while protecting caregiver safety.



TRANSFER




TURN



REPOSITION



PRONE



30° Positioning  
Wedges

HT-Air® 1200  
Air Supply

HoverCover™  
Breathable Absorbent  
Pad

The HoverMatt is  
X-ray friendly &  
radiolucent

The HoverMatt  
reduces the force  
needed to move a  
patient by 80–90%

HoverMatt® SPU Link

To reduce the need for boosting, the HoverMatt SPU Link comes equipped with unique connecting straps that attach to the bed to keep the device in place.





# HOVERSLING®

Repositioning Sheet

One Product, Many Uses

The HoverSling® Repositioning Sheet functions as both a HoverMatt SPU and sling to streamline patient handling tasks.

## AS A SLING



Sling Straps

True supine lift with reduced hip flexion

Sling Strap Pockets

## WITH AIR ASSIST



Transfer Handles (on both sides)

HT-Air® 1200 Air Supply



Adding Positioning Wedges provides the optimal solution for pressure injury prevention

- 30° positioning complies with NPIAP guidelines
- No-slip cover to secure and stabilize the device under the patient
- Available in disposable, reusable, and bariatric (reusable) options



SEATED LIFT/  
TRANSFER



SUPINE LIFT/  
TRANSFER



TURN



PRONE

# Q2ROLLER®

Staying within the Q2 protocol has never been easier & safer.



# HOVERJACK®

Lift fallen individuals from the floor to a bed or stretcher with no manual effort.



**Pro Tip:** Use a HoverMatt to transfer the fallen individual onto the HoverJack, then to the bed/stretcher.



Visit [www.HoverMatt.com](http://www.HoverMatt.com)  
for more information.



4482 Innovation Way  
Allentown, PA 18109

800.471.2776  
Fax 610.694.9601  
Info@HoverMatt.com

BRO\_ICU\_0321