



4482 Innovation Way
Allentown, PA 18109
800.471.2776
salesinfo@hovermatt.com
www HoverMatt.com

PATIENT SLING COMPATIBILITY STATEMENT

Information for Using HoverTech International HoverSling® Split-Leg/HoverSling® Repositioning Sheet with other Manufacturers' Lift Equipment.

HoverTech International has designed and manufactured safe patient handling products since its inception in 1997. HoverTech International is registered with the Food and Drug Administration and is subject to periodic audits by the FDA to verify our compliance with Good Manufacturing Practices. HoverTech International governs its own research, design and manufacturing following the guidelines and testing requirements set out in ISO 10535:2021 - Hoists for the transfer of persons: Requirements and Test Methods. HoverTech's slings are sold through national distributors and manufacturers of patient lifting equipment and are routinely used on lift equipment that is not manufactured by HoverTech International.

HoverSling® Split-Leg/HoverSling® Repositioning Sheets are fully compatible with other manufacturers' floor and ceiling lift systems using hook style hanger bars. Lifting slings may be used with 2, 3 or 4 point hook style hanger bars. Repositioning slings may be used with multiple point hook style hanger bars suitable for the size of sling being used. The care provider must follow HoverTech's User Manuals and recommendations with respect to the use, care, cleaning and inspection of slings and patient lifts.

COMPATIBLE

- Loop style Attachment
- Loop attachment on hook style-loop sling compatible hanger bars

NOT COMPATIBLE

- Clip style attachment Sling