



HOVERMATT®
PROS™ Air

HOVERMATT®
PROS™ Air Sling

Patient
Repositioning
Off-Loading
System

SYSTEM APPROACH TO PRESSURE INJURY PREVENTION

TARGETING:

- **Moisture:** HoverCover™; Breathable matt/sling.
- **Mobility:** Sling to transfer to chair/recliner;
Kitted with HoverMatt® PROSWedge™ to off-load pressure.
- **Friction & Shear:** Link straps to reduce matt/sling migration;
Air-assisted to reduce friction in movement.

HOVERMATT® PROS™ Air



PROS-HM-KIT shown includes PROS Air Matt with a HoverCover and pair of PROSWedges w/ & w/o tail.

HOVERMATT® PROS™ Air Sling



PROS-HMSL-KIT shown includes PROS Air Sling with a HoverCover and a pair of PROSWedges w/ & w/o tail.

BREATHABLE • COMPATIBLE • SINGLE PATIENT USE




Weight limit: 1200lbs

HOVERMATT® PROS™ Air

- Boosting
- Turning
- Repositioning
- Off-Loading w/ wedges
- Lateral Transfer
- Proning

PROS AIR & PROS AIR SLING COMMON FEATURES



- Kitted option with or without wedges



- Handles along both sides
- Extended turning handles



- Links to bed
- Full length
- Wipeable




- HoverCover holding patch
- QR Code instructions
- Head/foot indicators



- Folded for ease of placement



BREATHABLE • COMPATIBLE • SINGLE PATIENT USE



HOVERMATT® PROS™ Air Sling

- Boosting
- Turning
- Repositioning
- Proning
- Off-Loading w/ wedges
- Vertical Lifts:
True Supine & Seated
- Lateral Transfer

Weight limit:

- 600lbs as Lift
- 1200lbs as Air Assist

Product Numbers



HoverMatt® PROS™ Air

PROS-HM-KIT (Air Matt + HoverCover + pair of wedges*) 3 per case
PROS-HM-CS (Air Matt + HoverCover) 10 per case

*Wedge pair includes: 1 Wedge with Tail & 1 without Tail, Compressed



HoverMatt® PROS™ Air Sling

PROS-HMSL-KIT (Air Sling + HoverCover + pair of wedges*) 3 per case
PROS-HMSL-CS (Air Sling + HoverCover) 5 per case

*Wedge pair includes: 1 Wedge with Tail & 1 without Tail, Compressed



HoverMatt® PROSWedge™

PROS-WT (Pair of Wedges W/ & W/O Tail, Compressed) – 5 pairs per box
PROS-W (Pair of Wedges W/O Tail, Compressed) – 5 pairs per box



Full Size HoverCover™

HMAC-100-CASE (50 per case)



HoverMatt® PROS™ and HoverMatt® PROS™ Sling

Non-air slide sheet and sling w/ HoverCover also available kitted w/wo wedges.

